

Power Battery

WeOxy

Extended-runtime power battery designed for the WeOxy Q5, Q6 & Q7 concentrators, ideal for full-day mobility.

SPECIFICATIONS

Item Number	BA-704
Battery Type	Rechargeable Lithium-ion (Li-ion)
Voltage	14.4V
Rated Capacity	6.9Ah (99.3Wh)
Cell Model	41NR19/66-2
Certifications	CE, FCC
Runtime	Up to 9 hours on flow setting 1
Charge Time	Approx. 2 hours with included charger
Compatibility	WeOxy Q5, Q6 & Q7

WARRANTY

1 Year limited warranty. Covers manufacturing defects under normal use conditions.

This battery is intended for use only with designated WeOxy portable oxygen concentrators. Use of non-approved accessories may void the warranty. For replacement or support, contact support@weoxy.com.