

WeOxy Q7

Portable Oxygen Concentrator

WeOxy



PRODUCT SPECIFICATIONS

Model	WeOxy Q7
Dimensions	7 × 3 × 7.6 in
Weight	4.12 lbs
Direct Strap-Mount Design	Heavy Duty (No carry bag needed)
Inspiratory trigger	Ultra sensitive
Flow settings	1, 2, 3, 4, 5, 6, 7 (S)
Oxygen concentration	Up to 95.5%
Oxygen output	1,470 mL/min (setting S)
Battery life	Up to 9 hrs
Charging time	2 hrs
Warm-up time	2 minutes
AC power input	100–240 VAC, 100VA, 50–60 Hz
DC power input	13.0–16.8 VDC, 100W
Display	2.4" HD full-color screen
Assigned HCPCS code(s)	E1390+E1392

S (Summit Mode): High-output setting for elevated oxygen needs – ideal for demanding activities or high-altitude conditions.

BATTERY RUNTIME

Flow Setting	Flow	Power Battery
1	210 mL/min	Up to 9:00
2	420 mL/min	Up to 6:30
3	630 mL/min	Up to 3:50
4	840 mL/min	Up to 2:35
5	1050 mL/min	Up to 2:10
6	1260 mL/min	Up to 1:45
7 (s)	1470 mL/min	Up to 1:25

WHAT'S IN THE BOX

- ✓ Portable Oxygen Concentrator
- ✓ Carry Bag
- ✓ Shoulder Strap
- ✓ AC Power Cord & Adapter
- ✓ DC Power Car Charger
- ✓ Nasal Cannula
- ✓ User Manual
- ✓ Quick Start Guide
- ✓ Air Intake Filter Replacement Kit

WARRANTY COVERAGE

3yr

Device

2yr

Sieve bed

1yr

Battery



FDA Cleared



FAA Approved



Request a test unit.

✉ info@weoxy.com

☎ 1-(888)-999-2132