

Standard Battery

WeOxy

Lightweight standard battery delivering reliable runtime for the WeOxy Q5 portable oxygen concentrator.

SPECIFICATIONS

Item Number	BA-408
Battery Type	Rechargeable Lithium-ion (Li-ion)
Voltage	14.4V
Rated Capacity	6.5Ah (93.6Wh)
Cell Model	41NR19/66-2
Certifications	CE, FCC
Runtime	Up to 7 hours on flow setting 1
Charge Time	Approx. 2 hours with included charger
Compatibility	WeOxy Q5

WARRANTY

1 Year limited warranty. Covers manufacturing defects under normal use conditions.

This battery is intended for use only with designated WeOxy portable oxygen concentrators. Use of non-approved accessories may void the warranty. For replacement or support, contact support@weoxy.com.