

for WeOxy Portable Oxygen Concentrators

WeOxy Portable Oxygen Concentrators are medical devices intended to provide supplemental oxygen as prescribed by, or under the direction of, a licensed healthcare professional.

To support safe use during air travel, WeOxy Portable Oxygen Concentrators are designed to conform to applicable FAA acceptance criteria for portable oxygen concentrators intended for carriage and use on board aircraft.

On May 24, 2016, the Federal Aviation Administration (FAA) and the Department of Transportation (DOT) published amendments related to portable oxygen concentrators under 14 CFR 1.1, 1.2, 121.574, 125.219, and 135.91. The final rule, titled Oxygen and Portable Oxygen Concentrators for Medical Use by Passengers, replaced the previous model-by-model approval process and allows passengers to use a portable oxygen concentrator on board an aircraft if the device satisfies applicable FAA acceptance criteria and bears a label indicating conformance with those criteria.

FAA Compliance Statement Displayed on the Device



The manufacturer of this POC has determined this device conforms to all applicable FAA acceptance criteria for POC carriage and use on board aircraft.

FAA-Compliant Devices

Description	Model
WeOxy Portable Oxygen Concentrator	WeOxy Q5/WeOxy Q6/WeOxy Q7

Travel Advisory

- Before traveling, passengers should present the FAA compliance label on the device to airline personnel if requested.
- Passengers should contact the airline in advance to confirm any airline-specific requirements for the use of portable oxygen concentrators during travel. Requirements may vary by airline, and some airlines may request additional documentation or advance notification.
- **At least 150% of the expected maximum flight duration** may need to be covered by the batteries passengers bring to power their portable oxygen concentrator, as permitted under DOT requirements in 14 CFR Part 382. Passengers should consult their airline before travel to confirm specific battery and documentation requirements.
- Please follow all instructions in the applicable WeOxy User Manual, including battery charging, device operation, and safety precautions during travel.